



Greening Up The Lab: *10 Takeaway Tips*

1 Keep freezers clear of frost build-up and dispose of outdated samples.

2 Keep refrigerator/freezer filters and coils clean and free of obstructions.

3 Substitute hazardous cleaning and sterilizing substances with non-hazardous ones.

4 Follow proper waste disposal procedures.

5 Close variable speed fume hoods when not in use

6 Consider green alternatives when replacing or upgrading equipment.

7 Follow manufacturer recommendations when replacing HVAC and water filters.

8 Control water flow with solenoid valves, timers and automatic shut-off valves.

9 Set up inventory procedures to allow for effective resource management.

10 Share resources when possible to avoid unnecessary purchases.

